

# CEILING MEDICAL PENDANT



CONFIGURATION SHEET

2024




ARCHITECTURAL








ARIGMED



**SELECT THE NUMBER OF SECTIONS OF THE MEDICAL ARM:**



**SERVICE COLUMN SELECTION:**



**MEDICAL PENDANT SYSTEM QUANTITY: AREA DIMENTION:**

1  2  3 Medical Pendant quantity: \_\_\_\_\_

Height (ceiling - floor): \_\_\_\_\_ m.

Length room: \_\_\_\_\_ m.

**SPECIAL ELECTRICAL, DATA AND LIGHTING DEVICES:**

Ground Power Module  RJ45 Data inlet  Navigation Lighting  Service Column Lighting

**NURSE CALL PREPARATION:**

High: \_\_\_\_\_ cm

Width: \_\_\_\_\_ cm

**ELECTRICAL ARNESS:**

AWG 10  THW

AWG 12  XLP

HALOGEN FREE

**MEDICAL GAS OUTLETS:**

QUANTITY:

Oxygen Pieces: \_\_\_\_\_

Med Air Pieces: \_\_\_\_\_

Vacuum Pieces: \_\_\_\_\_

WAGD Pieces: \_\_\_\_\_

Nitrous Oxide Pieces: \_\_\_\_\_

Carbon Dioxide Pieces: \_\_\_\_\_

Nitrogen Pieces: \_\_\_\_\_

Instrumental Air Pieces: \_\_\_\_\_

**ELECTRICAL OUTLETS (HOSPITAL GRADE NEMA 5-20 , 110-127VAC 20 AMP):**

Red Pieces: \_\_\_\_\_  White Pieces: \_\_\_\_\_

Beige Pieces: \_\_\_\_\_  Ground Jack Pieces: \_\_\_\_\_

Orange ground fault Pieces: \_\_\_\_\_  Half turn 220v Pieces: \_\_\_\_\_

**MEDICAL GAS OUTLET CONNECTION TYPE (ARIGMED ´S ONLY):**

PURITAN BENNET  PLARRE type (only to Mexico)

MC70 AGA  CHEMETRON+OHMEDA (Diamond) AGA

DISS CGA-V5  MC70+M24 screw (only to Mexico)

WITHOUT (HOSES AND SUPPORTS ARE NOT INCLUDE)  CHEMETRON

**ACCESORIES TO INCLUDE:**

Medical basket Quantity: \_\_\_\_\_  Vacuum bottle support Quantity: \_\_\_\_\_  Infusion Pump support Quantity: \_\_\_\_\_  Tray Quantity: \_\_\_\_\_  Drawer Quantity: \_\_\_\_\_

I.V. hook support Quantity: \_\_\_\_\_  Vertical Rail Quantity: \_\_\_\_\_  2 I.V. + 2 Infusion Pump support system Quantity: \_\_\_\_\_  1 I.V. + 1 Infusion Pump support system Quantity: \_\_\_\_\_

**AVISO LEGAL PROPIEDAD INTELECTUAL E INDUSTRIAL:** Todos los derechos están reservados. Este documento es de la propiedad intelectual en contenidos e imágenes de **Grupo Arigmed S de R.L. de C.V.**, cualquier tipo de copia, similitud o plagio será objeto de las acciones legales pertinentes para evitar el uso por personas u organizaciones no autorizadas, así como resarcir el daño ocasionados al autor.



CERTIFICATE:

